

Invacare® 1.0L/1.7L Cylinder with Dual Lumen Conserver

For HomeFill® Oxygen System

en **HF2ECE6, HF2ECE9**
User Manual

This manual **MUST** be given to the user of the product.
BEFORE using this product, read this manual and save for future reference.



Yes, you can.®

©2017 Invacare Corporation

All rights reserved. Republication, duplication or modification in whole or in part is prohibited without prior written permission from Invacare. Trademarks are identified by ™ and ®. All trademarks are owned by or licensed to Invacare Corporation or its subsidiaries unless otherwise noted.

Making Life's Experiences Possible is a registered trademark in the U.S.A.

WD-40 and 3-IN-ONE are registered trademarks of the WD-40 Company.

Contents

1 General	4	6 After Use	25
1.1 Symbols	4	6.1 Disposal	25
1.2 Intended Use	5	6.2 Wear and Tear	25
1.3 User Profile	5	7 Troubleshooting	26
1.4 Service Life	5	7.1 Troubleshooting	26
1.5 Warranty	5	8 Technical data	27
1.6 Component Identification	6	8.1 Technical Description	27
2 Safety	7	8.2 Specifications	28
2.1 Label Locations	7		
2.2 General Guidelines	8		
3 Setup	14		
3.1 Dealer Setup	14		
3.2 User Setup	14		
3.3 Receiving Inspection	14		
3.4 Cylinder Prefill Inspection	14		
3.4.1 The Hydrostatic Testing Date	15		
3.4.2 External Examination	16		
4 Usage	17		
4.1 Before Use	17		
4.2 Connecting/Disconnecting the Cylinder to/from the Compressor	17		
4.3 Cylinder Assembly Operation	17		
5 Maintenance	21		
5.1 Maintenance	21		
5.2 Emptying the Cylinder Assembly	22		
5.3 Cleaning	23		
5.4 Disinfecting	23		
5.5 Reuse	23		
5.6 Storage	24		
5.7 Optional Accessories	24		

1 General

1.1 Symbols

Signal words are used in this manual and apply to hazards or unsafe practices which could result in personal injury or property damage. See the information below for definitions of the signal words.



DANGER!

– Danger indicates an imminently hazardous situation which, if not avoided, will result in death or serious injury.



WARNING!

– Warning indicates a potentially hazardous situation which, if not avoided, could result in death or serious injury.



CAUTION!

– Caution indicates a potentially hazardous situation which, if not avoided, may result in property damage or minor injury or both.



– Indicates a hazardous situation that could result in damage to property if it is not avoided.



Gives useful tips, recommendations and information for efficient, trouble-free use.



Read the manual



No open flame



DO NOT Smoke



Oxidizer



Recycle



DO NOT dispose of in household waste



Manufacturer



EU Representative



This product complies with Directive 93/42/EEC concerning medical devices.

The launch date of this product is stated in the CE declaration of conformity.

1.2 Intended Use

This cylinder is intended to store, regulate, and deliver supplemental oxygen to patients who may have difficulty extracting oxygen from the air they breathe. It is for patients who would normally receive the oxygen via a nasal cannula. The device delivers USP 93±3% oxygen at a number of different flow settings and may increase the drying out of airways. This product is to be used as an oxygen supplement and is not intended to be life supporting or life sustaining.



DANGER! **Risk of Injury or Death**

This product is to be used as an oxygen supplement and is not intended to be life-supporting or life-sustaining. ONLY use this product if the patient is capable of spontaneous breath, able to inhale and exhale without the use of a machine.



WARNING! **Risk of Injury or Damage**

Use of this product outside of the intended use and specifications has not been tested and may lead to product damage, loss of product function, or personal injury.

- DO NOT use this product in any way other than described in the specifications and intended use sections of this manual.

1.3 User Profile

No specific product knowledge or training is required, other than what is contained in this manual, to operate the

1176148-D

product. It is not intended for supplying oxygen therapy during refilling. This product is intended to be used by or under adult supervision only after reading and understanding the instructions and warnings of this user manual.

Only personnel instructed and trained in the use of the HomeFill® Conserving Cylinder should operate this product.

The use of this product is limited to the oxygen patient. Cylinders MUST be used only by the oxygen patient and are not to be distributed to any other individual for any purpose.

This product is NOT to be used by patients who breathe through their mouths.

1.4 Service Life

The expected service life of the conserver is five years when used in accordance with the safety instructions, maintenance intervals and correct use stated in this manual. Refer to the Maintenance section.


The expected service life of the cylinder is based on testing the cylinder every ten years and using it in accordance with the safety instructions, maintenance intervals and correct use stated in this manual. Refer to the following procedures:

- The Hydrostatic Testing Date, External Examination, Maintenance

1.5 Warranty

Terms and conditions of the warranty are part of the general terms and conditions particular to the individual countries in which this product is sold.

1.6 Component Identification



Item	Description
A	Cylinder Assembly
B	Cylinder
C	Continuous Flow Conservor
D	Flow Dial
E	Indicator Window
F	Switch
G	Outlet Connection
H	Cylinder Fillport
I	Shutoff Valve
J	Cylinder Gauge
K	Cylinder Gauge Needle

2 Safety

2.1 Label Locations




Table of regulatory information with columns for country codes (GB, CZ, DK, D, E, H, F, I, NL, N, PL, S) and corresponding safety text in English and other languages.

i Serial number located above product label shown here.

2.2 General Guidelines

In order to ensure the safe installation, assembly and operation of the product these instructions MUST be followed.



DANGER!

Risk of Death, Injury or Damage

- Improper use of this product may cause death, injury or damage. This section contains important information for the safe operation and use of this product.
- DO NOT use this product or any available optional equipment without first completely reading and understanding these instructions and any additional instructional material such as user manuals, service manuals or instruction sheets supplied with this product or optional equipment.
- If you are unable to understand the warnings, cautions or instructions, contact a healthcare professional, dealer or technical personnel before attempting to use this equipment.
- Check ALL external components and carton for damage. In case of damage, or if the product is not working correctly, contact a technician or Invacare for repair.
- THE INFORMATION IN THIS DOCUMENT IS SUBJECT TO CHANGE WITHOUT NOTICE.



DANGER!

Risk of Death, Injury or Damage

- Use this cylinder in accordance with the HomeFill compressor user manual.
- DO NOT handle the cylinder assembly or use the contents until you are professionally trained, including emergency procedures. Training required to operate the product is contained in the user manual.
- This product is not intended for supplying oxygen therapy during refilling.
- This product is intended to be used by or under adult supervision.
- Oxygen Cylinders filled by the Invacare HomeFill oxygen system shall be used for personal use only. “Not to be filled for resale or use by professional users.”
- This cylinder is to be used with the Invacare HomeFill oxygen system. This cylinder contains 93±3% oxygen produced by an oxygen concentrator using the pressure swing adsorption process.
- ALWAYS open the valve slowly and close the valve after each use and when empty.

**DANGER!****Risk Of Death, Injury, Or Damage From Fire**

Textiles, oil or petroleum substances, grease, greasy substances and other combustibles are easily ignited and burn with great intensity in oxygen enriched air and when in contact with oxygen under pressure. To avoid fire, death, injury or damage:

- DO NOT SMOKE while using this product.
- DO NOT use near OPEN FLAME or IGNITION SOURCES.
- DO NOT use any lubricants on the cylinder unless recommended by Invacare.
- NO SMOKING signs should be prominently displayed.
- Avoid creation of any spark near oxygen equipment. This includes sparks from static electricity created by any type of friction.
- Keep all matches, lighted cigarettes, electronic cigarettes or other sources of ignition out of the room in which this cylinder is located and away from where oxygen is being delivered.
- Keep the oxygen tubing, cord, and cylinder out from under such items as blankets, bed coverings, chair cushions, clothing, and away from heated or hot surfaces including space heaters, stoves, and similar electrical appliances.

**CAUTION!**

Federal (statutory) law restricts this device sale to or on the order of a medical practitioner licensed by a governmental agency where he/she practices.

- ONLY a licensed medical practitioner may order the purchase or use of this device.

Operating Information



DANGER! **Risk of Fire**

Dirt or contaminants on the cylinder assembly or conservor can cause fire, injury or death.

- DO NOT use if dirt or contaminants are present on or around cylinder or conservor.
- DO NOT use or fill cylinder assembly if the cylinder fillport is dirty or contaminated. Contact your provider for service or repair.
- DO NOT allow oil from your hands or other sources to come into contact with the compressor coupling or cylinder assembly fittings. These substances may become flammable in the presence of oxygen and cause injury. Examples of items to avoid are hand cream or lotion, cooking oil, suntan oil or sunblock and similar oily products. If you use these types of products Invacare suggests you wash your hands prior to using the HomeFill system. If the cylinder fillport should come in contact with such substances, contact your provider for service. Avoid touching the cylinder fillport.



DANGER! **Risk of Fire**

Dirt or contaminants on the cylinder assembly or conservor can cause fire, injury or death.

- DO NOT clean with aromatic hydrocarbons, including, but not limited to benzene and toluene.
- NEVER oil or lubricate the compressor coupler or cylinder fillport. NEVER use any penetrating oil or lubricant such as WD-40® or 3-IN-ONE®. If connection difficulty is experienced, verify that the cylinder fillport is being inserted straight down into the compressor coupler. If continued difficulty is experienced, contact your provider for assistance.



DANGER!
Risk of Fire

Dirt or contaminants on the cylinder assembly or conserver can cause fire, injury or death.

- DO NOT use cleaning/disinfecting solutions on the fillport fitting.
- DO NOT immerse product in any kind of liquid.
- DO NOT allow any kind of liquid to get inside openings of the conserver.
- DO NOT spray cleaning/disinfecting solutions on product.
- DO NOT autoclave.
- DO NOT gas sterilize with ethylene oxide.
- DO NOT attempt to repair cylinder assembly. All repairs MUST be performed by Invacare Corporation or in an Invacare approved repair center.
- Store cylinder assembly in a clean, dry and cool area when not in use.



WARNING!
Risk of Injury or Death

Contaminants in the oxygen may cause injury or death.

- DO NOT use if the oxygen has an unusual odor.
- DO NOT use if particulates are inhaled.
- Contact Invacare or a dealer for service.



WARNING!
Risk of Injury or Damage

To avoid injury or damage from product misuse:

- The flow dial should be turned with normal pressure from one's hand. If the flow dial is stiff to turn, do not force, otherwise damage or injury could occur.
- DO NOT use a flow outlet for driving any medical equipment.



DANGER!
Risk Injury or Death

To avoid choking or ingestion of chemicals from airway contamination:

- DO NOT use the product in the presence of pollutants, smoke, fumes, flammable anesthetics, cleaning agents, or chemical vapors.



CAUTION!

Risk of Injury or Damage

To reduce the risk of injury or damage from dropped or impacted cylinder assemblies:

- Use two hands when handling/transporting oxygen cylinders.
- Only use Invacare Corporation carrying bag.
- DO NOT allow cylinder assemblies to tip or fall. DO NOT drop cylinder assemblies.
- Secure gas cylinders so they cannot fall.
- Ensure cylinder is always in upright position when not connected to the HomeFill compressor.
- DO NOT allow cylinders to impact inanimate objects.



CAUTION!

Risk of Damage

To avoid product damage:

- DO NOT block the outlet fitting or kink the cannula when the product is in use, as this may damage the unit.
- DO NOT fill the cylinder assembly outdoors or anywhere below 13°C (55°F).



DANGER!

Risk of Injury or Death

The cannula can cause tripping, falling, or other injury if improperly positioned and secured. To avoid injury or death:

- The cannula MUST be routed and secured properly.
- DO NOT position the cannula around the neck. Ensure the patient can move freely while wearing the cannula.
- Avoid positioning cannulas across areas of high foot traffic (i.e. aisles, doorways, hallways, etc.)



WARNING!

Risk of Injury or Death

To avoid choking and/or strangulation from tubing entanglement:

- Keep children and pets away from nasal cannula and tubing.
- Close supervision is necessary when the nasal cannula is used by or near children and/or disabled persons.

**WARNING!****Risk of Injury or Death**

To reduce the risk of injury or death from illness:

- Before the first patient use, the provider **MUST** disinfect the product. Refer to the disinfecting section.
- Replace the nasal cannula on a regular basis. Check with your equipment provider or physician to determine how often the cannula should be replaced.
- **DO NOT** share cannulas between patients.



The cylinder poses no biocompatibility hazards and is latex and phthalate free. The biocompatibility of cylinder accessories may vary by manufacturer. Invacare accessories such as the cylinder bag and Invacare oxygen tubing are latex and phthalate free. Refer to the packaging on any cylinder accessories to determine the biocompatibility hazard.

The cylinder assembly **MUST** be emptied of gas prior to shipping by courier or postal service. Refer to Emptying the Cylinder Assembly.

Because the cylinder contains magnetic, ferrous material, it is not suitable to be in an MRI room.

Be sure all connections are tight and leak free. **DO NOT** use a liquid leak detector to test for leaks. Leaks are indicated by loss of pressure when the cylinder assembly is not in use and by hissing sounds. This product is equipped with a relief valve. If you hear a loud hissing or popping sound coming from the product, discontinue use of product, turn the flow dial to the “0” setting, and contact your dealer/Home Care equipment provider.

ALWAYS confirm prescribed dose before administering to patient and monitor on a frequent basis.

NO OXYGEN is delivered when the conserver flow dial setting is “0”.

DO NOT use a humidifier with this product.

Keep cylinder in a well ventilated place.

This product is designed to operate with a standard or high flow dual lumen, adult oxygen nasal cannula with a maximum length of 2.1 m (7 ft).

The cannula is for single patient use only.

Cannula **MUST** be installed on outlet connection prior to use.

Refer to the HomeFill User Manual for additional warnings, cautions and safety instructions.

3 Setup

3.1 Dealer Setup



WARNING!

- Before the first patient use, the provider **MUST** disinfect the product. Refer to the Disinfecting section.
- Before patient use, the provider **MUST** check the hydrostatic testing date. Refer to the Hydrostatic Testing Date section.



WARNING!

Risk of Injury

- For medical applications, it is very important to select the prescribed level of oxygen flow. Uninterrupted use of high concentrations of oxygen over a long duration, without monitoring its effect on oxygen content of arterial blood, may be harmful.
- Use **ONLY** as directed by a licensed practitioner.
 - **DO NOT** increase or decrease the flow unless a change has been prescribed by your physician or therapist.

3.2 User Setup



WARNING!

- When receiving the cylinder from the provider, the user/patient **MUST** perform a receiving inspection. Refer to the Receiving Inspection section.

3.3 Receiving Inspection

Remove product from package and inspect for damage. If there is any damage, **DO NOT** use. Contact your dealer/Home Care equipment provider.

3.4 Cylinder Prefill Inspection



WARNING!

Risk of Injury or Damage

- Oxygen cylinders shall be used for personal use only. Not to be filled for resale or use by professional users.
- Filling of cylinders is performed by Invacare HomeFill oxygen system only.
- This cylinder is to be used with the HomeFill system. This cylinder contains 93±3% oxygen produced by an oxygen concentrator using the pressure swing adsorption process.
- **DO NOT** handle the cylinder or use contents until you are professionally trained, including emergency procedures.
- Use this cylinder in accordance with the HomeFill user manual.
- Open and close the shutoff valve slowly. Close the shutoff valve after each use and when empty.
- Secure the cylinder during storage and use. **DO NOT** drop.
- No smoking in cylinder area. Keep away from heat, flame or spark.
- Keep out of the reach of children.

**DANGER!****Risk of Injury or Death**

Improper filling or abuse of this cylinder, or failure to follow these safety instructions may cause serious bodily injury or death. High pressure oxidizing gas vigorously accelerates combustion.

- DO NOT alter or modify cylinder or related components.
- DO NOT use oil or grease. Keep cylinder, shutoff valve, all equipment and connections free of oil and grease to avoid violent ignition resulting in serious bodily injury or death.

3.4.1 The Hydrostatic Testing Date

Aluminum cylinders MUST undergo testing every ten years.

**DANGER!****Risk of Death, Injury or Damage from Fire**

To avoid fire, death, injury or damage:

- Cylinders that have been refinished and/or exposed to elevated temperatures MUST be hydrostatically tested before filling. Cylinder temperature should not exceed 52° C (125° F). Aluminum cylinders subject to the action of fire or heated temperatures in excess of 177° C (350° F) MUST be withdrawn from service and condemned by trained personnel. DO NOT use caustic paint strippers. Valve and safety relief devices MUST be removed and replaced by trained and authorized personnel. DO NOT alter or change valve assembly.
- All cylinders MUST be inspected before attempting to fill; otherwise, injury or damage may occur.

**WARNING!****Risk of Injury**

To reduce the risk of injury:

- DO NOT fill cylinders that have not been tested in the past ten years. Otherwise, serious personal injury may result. Contact your dealer/Home Care equipment provider for replacement.

3.4.2 External Examination



WARNING!

Risk of Injury or Damage

To reduce the risk of injury or product damage:
– DO NOT use the cylinder if any of the following items are found. Contact your dealer/Home Care equipment provider for replacement.

1. Examine the outside of the cylinder for the following conditions. If they exist, contact your provider for service.
 - Dents or dings, arc burns, oil or grease, any other signs of damage that might cause a cylinder to be unacceptable or unsafe for use
2. Examine the cylinder for evidence of fire or thermal damage. Evidence includes charring or blistering of the paint, or other protective coating or heat sensitive indicator. If fire or thermal damage is found, contact your dealer/Home Care equipment provider for a replacement.
3. Inspect the cylinder fillport for the following. If they exist, contact your provider for a replacement.
 - Debris, oil or grease, noticeable signs of damage, signs of corrosion inside the valve, signs of excessive heat or fire damage, signs of pressure leaks

4 Usage

4.1 Before Use



WARNING!

Risk of Injury or Damage

Defective or inoperable cylinders may cause injury or damage.

- All cylinders **MUST** be inspected before attempting to fill.
- If the cylinder is found to be defective or does not work properly, contact an Invacare dealer or qualified technician for service or replacement.

Before filling a cylinder, read and follow the information in the following sections:

- Cylinder Prefill Inspection
- External Examination

4.2 Connecting/Disconnecting the Cylinder to/from the Compressor



Refer to the HomeFill Compressor user manual for operating instructions to connect and disconnect the cylinder to the HomeFill Compressor.


4.3 Cylinder Assembly Operation



WARNING!


Risk of Fire

- NEVER smoke in an area where oxygen is being administered.
- NEVER use near any type of flame or flammable/explosive substances, vapors or atmosphere.



Invacare® 1.0L/1.7L Cylinder with Dual Lumen Conservor

1. Prior to each use, inspect product for visible damage. Refer to the External Examination section. DO NOT use if any damage is found.

 If for any reason, any label becomes illegible or lost, contact your dealer/Home Care equipment provider.

2. Inspect the cylinder gauge ① to ensure the pressure is at least 1500 psig (10342 kPa). If below 1500 psig (10342 kPa), then refill the cylinder ② using the HomeFill system.



WARNING!

– When cylinder gauge needle on cylinder gauge enters red zone, it is recommended to refill the cylinder on the Invacare HomeFill Compressor.

3. Position the cylinder so the oxygen cylinder valve outlet ③ is pointing away from the user and any other person(s).



WARNING!
Risk of Injury

To reduce the risk of injury:
– DO NOT direct flow of oxygen at any person, or flammable material when adjusting the flow dial.


4. Turn the flow dial ④ to the “0” setting.

5. Attach a dual lumen adult oxygen nasal cannula with a maximum length of 2.1 m (7 feet) to the product’s two outlet connections according to the cannula manufacturer’s instructions.



It does not matter which tube is installed on which outlet connection.

DO NOT use pediatric, low flow nasal cannulas or oxygen masks with this product.




DANGER!

Risk of Death

The cannula can cause tripping, falling, or other injury if improperly positioned and secured. To avoid injury or death:


- The cannula MUST be routed and secured properly.
- DO NOT position the cannula around the neck. Ensure the patient can move freely while wearing the cannula.
- Avoid positioning cannulas across areas of high traffic (i.e. aisles, doorways, hallways, etc.)


6. Place the cannula over your ears and position the prongs in your nose as instructed by your Health Care provider or cannula manufacturer.
7. Ensure there are no kinks or twists in the cannula tubing.
8. Open the Shutoff valve ① knob slowly towards the ON direction (counterclockwise)

 This will allow cylinder pressure to travel through the converter.


9. Use the switch ② to select one of the following modes:

- Continuous Flow ③
- Conserving Pulse ④

 The converter is labeled “CONTINUOUS” and “PULSE”. The first three flows (0.25, 0.50 and 0.75 L/min) are in continuous mode at all times, regardless of the switch position.



10. Turn flow dial until the indicator window ⑤ is aligned with the prescribed setting.

 When there are two numbers displayed in the indicator window, the larger number is the PULSE flow setting and the smaller number is the CONTINUOUS flow setting.



WARNING!

- To avoid injury to patient, ALWAYS confirm prescribed setting before administering to patient and monitor flow on a frequent basis.
- Only use the Invacare Corporation carrying bag.

Invacare® 1.0L/1.7L Cylinder with Dual Lumen Conserver

11. Breathe through the nose and feel a pulse of oxygen as you inhale.



When on a pulse setting, there is flow, or a pulse, only at the beginning of each breath. If you do not feel the pulse at the beginning of each breath, check the setting.

Because each patient's breathing pattern is different and the environment varies, it may be difficult to feel some low setting pulses.

If operating the cylinder below 1.7 ° C (35° F) use the 2 L/Min continuous setting only.

When the unit is first switched from pulse to continuous flow, a single inhalation is required to initiate the continuous flow.

12. When finished using the cylinder, perform the following steps:
 - a. Turn the flow dial to the "0" setting.
 - b. Slowly close the shutoff valve (clockwise) fully.

5 Maintenance

5.1 Maintenance



WARNING!

Risk of Injury or Damage

Lack of product maintenance can lead to injury and loss of product function or damage.

- Maintenance **MUST** be performed **ONLY** by qualified personnel.
- Regular maintenance is necessary to assure proper operation.



WARNING!

Risk of Injury or Damage

– The cylinder **MUST** be disconnected from the HomeFill and emptied of gas prior to cleaning or disinfecting. Refer to the Emptying the Cylinder Assembly section. The fillport cover **MUST** also be installed.

– All cylinders **MUST** be inspected before attempting to fill; otherwise, injury or damage may occur. Refer to the Cylinder Prefill Inspection section.



DANGER!

Risk of Fire

Dirt or contaminants on the cylinder or conserver can cause fire.

- **DO NOT** use or fill cylinder if the cylinder fillport is dirty or contaminated. Contact your provider for service or repair.
- **DO NOT** clean with aromatic hydrocarbons, including, but not limited to benzene and toluene.
- **DO NOT** use cleaning/disinfecting solutions on the fillport fitting.
- **DO NOT** immerse product in any kind of liquid.
- **DO NOT** allow any kind of liquid to get inside openings of the conserver.
- **DO NOT** spray cleaning/disinfecting solutions on product.
- All repairs **MUST** be performed by Invacare Corporation or in an Invacare approved repair center.

The cannula does not require maintenance and needs to be replaced at the manufacturer's recommended interval.

The aluminum cylinder must undergo hydrostatic testing every 10 years. Refer to The Hydrostatic Testing Date section.

5.2 Emptying the Cylinder Assembly



DANGER!

Risk of Fire


- NEVER smoke in an area where oxygen is being administered.
- NEVER use near any type of flame or flammable/explosive substances, vapors or atmosphere.
- DO NOT direct flow of oxygen at any person, or flammable material when adjusting the flow dial.







CAUTION!

Risk of Damage


- DO NOT completely empty the cylinder. Otherwise, damage to the cylinder could occur.



1. Position the cylinder assembly **A** so the outlet connection **G** is pointing away from the user and any other person(s).
2. Turn the flow dial **D** to setting “0”.
3. Slowly open the shutoff valve **I** fully.
 - i** Turn the shutoff valve **I** counterclockwise to open.
4. Turn the flow dial **D** to setting “6”.
5. Examine the cylinder gauge **J** until the pressure is below 500 psig (3447 kPa).
 - i** Cylinder pressure is below 500 psig (3447 kPa) when the cylinder gauge needle **K** is in the red zone.

6. Turn the flow dial  to setting "0".
7. Close the shutoff valve  fully.
 -  Turn the shutoff valve clockwise  to close.

5.3 Cleaning


-  Refer to the safety information in Maintenance before performing this procedure.

Follow the cleaning/disinfectant product manufacturer's instructions for use, handling, storage and disposal of the product.

Follow this procedure to clean the cylinder as required when visibly soiled.

1. Ensure the cover is installed on the fillport.
2. Clean exterior surfaces of the cylinder with a lint free cloth dampened with 70% isopropyl alcohol.
3. Clean the exterior surfaces of the conserver with a dry, lint free cloth.
4. If the conserver is still visibly soiled, clean with a lint free cloth dampened with 70% isopropyl alcohol. DO NOT wipe the fillport or outlet fitting.
5. Wipe the cylinder and conserver with a dry lint free cloth. Ensure the product is dry before using.

5.4 Disinfecting

-  Refer to the safety information in Maintenance before performing this procedure.

Follow the cleaning/disinfectant product manufacturer's instructions for use, handling, storage and disposal of the product.

This procedure is for disinfecting the cylinder and is to be performed by the provider prior to a patient change and as otherwise required.

Recommended disinfectant products:

- Quaternary ammonium based wipes selected by the provider for the desired disinfecting characteristics
 - 70% isopropyl alcohol wipes
1. Ensure the cover is on the fillport.
 2. If heavily soiled, first clean the cylinder. Refer to Cleaning.
 3. Use a wipe to disinfect the external surfaces of the conserver and cylinder. Follow the disinfectant manufacturer's instructions. DO NOT wipe the fillport or outlet connection.
 4. Wipe the cylinder and conserver with a dry lint free cloth. Ensure the product is dry before using.

5.5 Reuse




WARNING!

Before transferring the product to a different user, providers MUST inspect the product for damage and disinfect the product. Refer to the following procedures:

- Disinfecting, Cylinder Prefill Inspection, The Hydrostatic Testing Date, External Examination

The cylinder can be transferred to a different patient. The cannula is for single patient use only.

5.6 Storage



1. Ensure the shutoff valve ① is closed.
2. Ensure the flow dial ② is turned to setting “0”.
3. Place the cover on the fillport ③.
4. Store product in a clean, dry and cool area free from grease, oil, and other sources of contamination.

5.7 Optional Accessories



WARNING!

Risk of Serious Injury or Damage

Use of non-Invacare accessories may result in serious injury or damage.

- Invacare products are specifically designed and manufactured for use in conjunction with Invacare accessories. Accessories designed by other manufacturers have not been tested by Invacare are not recommended for use with Invacare products.
- DO NOT use non-Invacare accessories.
- To obtain Invacare accessories, contact Invacare by phone or at www.invacare.com.

The following optional accessories are available for use with this product: HF2POST6BAG, HF2POST9BAG – Carrying Bag

6 After Use

6.1 Disposal

This product has been supplied from an environmentally aware manufacturer that complies with the Waste Electrical and Electronic Equipment (WEEE) Directive 2002/96/CE. This product may contain substances that could be harmful to the environment if disposed of in places (landfills) that are not appropriate according to legislation.

Follow local governing ordinances and recycling plans regarding disposal of the product or components normally used in operation. The device does not generate waste or residue in operation.

- DO NOT dispose of the device in the normal waste stream.
- Any accessories not part of the device MUST be handled in accordance with the individual product marking for disposal.

6.2 Wear and Tear

Invacare reserves the right to ask for any item back that has an alleged defect in workmanship. See Warranty that shipped with the product for specific warranty information.

Refer to the Maintenance section for proper preventative maintenance schedule and use of the product.

7 Troubleshooting

7.1 Troubleshooting

If Conserving Cylinder fails to function, consult the Troubleshooting guide. If problem cannot be corrected, consult your dealer/Home Care equipment provider.

Problem	Possible Cause	Solution
No flow.	<ol style="list-style-type: none"> 1. Flow dial at setting "0". 2. Cylinder empty. 3. Flow dial is between settings. 4. Kinked or blocked cannula. 5. Shutoff valve closed. 	<ol style="list-style-type: none"> 1. Adjust flow dial to prescribed setting. 2. Fill or replace cylinder. 3. Set dial so indicator points to a setting. 4. Check to see if the cannula is kinked. 5. Slowly open shutoff valve.
Cylinder will not fill.	<ol style="list-style-type: none"> 1. Shutoff valve closed. 	<ol style="list-style-type: none"> 1. Slowly open shutoff valve.

8 Technical data

8.1 Technical Description

The Invacare HF2ECE6 and HF2ECE9 Oxygen Conserving Cylinders are oxygen storage products that can deliver either continuous or pulse flow oxygen. The cylinders are designed to store compressed oxygen up to 14479 kPa (2100 psig). The continuous flow conserver reduces the pressure at the cannula connection to 345 kPa (50 psig) and contains a rotary flow dial that is used to adjust the oxygen flow rate to the prescribed setting. The cylinders are filled with USP 93±3% oxygen using their proprietary fillports that allow filling only through the use of an Invacare HomeFill System. The duration of oxygen flow between fills is dependent on the size of the cylinder and the flow setting.

8.2 Specifications

Pulse Setting	1, 1.5, 2, 2.5, 3, 4, 5, 6				
	Conserving (Savings) Ratio, Bolus Size and Accuracy at 20 BPM		Pulse Volumes (L)		
	Equivalent Flow Rate		Minimum	Nominal	Maximum
	1 L/min setting:		0.014	0.0165	0.019
	1.5 L/min setting:		0.021	0.0248	0.029
	2 L/min setting:		0.029	0.0330	0.038
	2.5 L/min setting:		0.036	0.0413	0.047
	3 L/min setting:		0.043	0.0495	0.057
	4 L/min setting:		0.056	0.0660	0.076
	5 L/min setting:		0.071	0.0825	0.095
	6 L/min setting:		0.084	0.0990	0.114
	Sensitivity		0.130 to 0.150 in (0.330 to 0.381 cm) H ₂ O		
	Delivery		"ON" time	less than 125 ms	
			"OFF" time	less than 300 ms	

Continuous Flow Setting	Pulse Volumes (L)		
	Equivalent Flow Rate	Minimum	Nominal
1/4 L/min setting:	0.213	0.250	0.288
1/2 L/min setting:	0.425	0.500	0.575
3/4 L/min setting:	0.638	0.750	0.863
0.4 L/min setting:	0.340	0.400	0.460
0.5 L/min setting:	0.425	0.500	0.575
1 L/min setting:	0.850	1.000	1.150
1.5 L/min setting:	1.275	1.500	1.725
2 L/min setting:	1.700	2.000	2.300
2.5 L/min setting:	2.125	2.500	2.875
3.5 L/min setting:	2.975	3.500	4.025
4.5 L/min setting:	3.825	4.500	5.175

Pressures

Inlet Pressure Range	2068–14479 kPa (300-2100 psig)
Max Fill Pressure	14479 kPa (2100 psig)
Internal Regulated Pressure	345 kPa (50 psig)
Safety Burst Disk	Bursts at 20.8–23.2 Mpa at 22°C (3025–3360 psig at 72°F)

Invacare® 1.0L/1.7L Cylinder with Dual Lumen Conserver

Environmental Operating Temperature Range

Storage Temperature Range	-20°C to 60°C (-4°F to 140°F)
Operating Temperature Range	1.7°C to 40.6°C (35°F to 105°F)
Max Operating Humidity	95% non-condensing


Assembly Dimensions

Model	HF2ECE6	HF2ECE9
Diameter	11.1 cm (4.4 in)	11.1 cm (4.4 in)
Weight	2.07 kg (4.57 lb)	2.75 kg (6.07 lb)
Overall Length	43.3 cm (17.03 in)	40.7 cm (16.03 in)

Cylinders

Cylinder Size	M6	M9
Volume (H₂O)	1.0 L (0.035 ft ³)	1.7 L (0.061 ft ³)
Oxygen Capacity	170 L (6 ft ³)	255 L (9 ft ³)

Fill Times

Model	HF2ECE6	HF2ECE9
	1 hr 25 min*	2 h 20 min*
* Concentrator flow rates to patient:		
Platinum 5, Platinum S Series and Perfecto ₂ Series	Up to 2.5 L/min	
Platinum 9	Up to 5 L/min	
 All filling times are approximate and may vary on environmental conditions.		

Accessories

Model	HF2ECE6	HF2ECE9
Dual Lumen Cannula	5' (1.524 m) or 7' (2.1336 m) tubing length (HUD 1850)	5' (1.524 m) or 7' (2.1336 m) tubing length maximum (HUD 1850)
Carry Bag	HF2POST6BAG	HF2POST9BAG

Regulatory Compliance

93/42/EEC - Medical Device Directive (Class 2b)
2014/68/EU - Pressure Equipment Directive (Category 3)

Notes

Notes

Notes

Belgium & Luxemburg: Invacare nv, Autobaan 22, B-8210 Loppem • Tel: (32) (0) 50 83 10 10 • Fax: (32) (0) 50 83 10 11 • belgium@invacare.com • www.invacare.be

Danmark: Invacare A/S, Sdr. Ringvej 37, DK-2605 Brøndby • Tel: (45) (0)36 90 00 00 • Fax: (45) (0)36 90 00 01 • denmark@invacare.com • www.invacare.dk

Deutschland: Invacare GmbH, Alemannenstraße 10, D-88316 Isny • Tel: (49) (0)75 62 7 00 0 • Fax: (49) (0)75 62 7 00 66 • kontakt@invacare.com • www.invacare.de

Eastern Europe, Middle East & CIS: Invacare EU Export • Kleiststraße 49 • D-32457 Porta Westfalica • Germany • Tel: (49) 5731 754540 • Fax: (49) 5731 754541 • webinfo-eu-export@invacare.com • www.invacare-eu-export.com

España: Invacare SA, c/Areny s/n, Polígon Industrial de Celrà, E-17460 Celrà (Girona) • Tel: (34) (0)972 49 32 00 • Fax: (34) (0)972 49 32 20 • contactsp@invacare.com • www.invacare.es

France: Invacare Poirier SAS, Route de St Roch, F-37230 Fondettes • Tel: (33) (0)2 47 62 64 66 • Fax: (33) (0)2 47 42 12 24 • contactfr@invacare.com • www.invacare.fr

Ireland: Invacare Ireland Ltd, Unit 5 Seatown Business Campus • Seatown Road, Swords, County Dublin – Ireland • Tel: (353) 1 810 7084 • Fax: (353) 1 810 7085 • ireland@invacare.com • www.invacare.ie

Italia: Invacare Mecc San s.r.l., Via dei Pini 62, I-36016 Thiene (VI) • Tel: (39) 0445 38 00 59 • Fax: (39) 0445 38 00 34 • italia@invacare.com • www.invacare.it

Nederland: Invacare AE, Galvanistraat 14–3, NL–6716 BZ Ede • Tel: (31) (0)318 695 757 • Fax: (31) (0)318 695 758 • nederland@invacare.com • csede@invacare.com • www.invacare.nl

Norge: Invacare AS, Grensesvingen 9, Postboks 6230, Etterstad, N-0603 Oslo • Tel: (47) (0)22 57 95 00 • Fax: (47) (0)22 57 95 01 • norway@invacare.com • island@invacare.com • www.invacare.no

Österreich: Invacare Austria GmbH, Herzog Odilostrasse 101, A-5310 Mondsee • Tel: (43) 6232 5535 0 • Fax: (43) 6232 5535 4 • info@invacare-austria.com • www.invacare.at

Portugal: Invacare Lda • Rua Estrada Velha, 949, P-4465-784 Leça do Balio • Tel: (351) (0)225 1059 46/47 • Fax: (351) (0)225 1057 39 • portugal@invacare.com • www.invacare.pt

Sverige: Invacare AB • Fagerstagatan 9 • S-163 53 Spånga • Tel: (46) (0)8 761 70 90 • Fax: (46) (0)8 761 81 08 • sweden@invacare.com • www.invacare.se

Suomi: Camp Mobility • Patamäenkatu 5, 33900 Tampere • Tel: 09-350 76 310 • info@campmobility.fi • www.campmobility.fi

Schweiz/Suisse/Svizzera: Invacare AG • Benkenstrasse 260 • CH-4108 Witterswil • Tel.: (41) (0)61 487 70 80 • Fax.: (41) (0)61 487 70 81 • switzerland@invacare.com • www.invacare.ch

United Kingdom: Invacare Limited, Pencoed Technology Park, Pencoed, Bridgend CF35 5HZ • Tel: (44) (0) 1656 776222 • Fax: (44) (0) 1656 776220 • UK@invacare.com • www.invacare.co.uk

Invacare Corporation

USA

One Invacare Way
Elyria, Ohio USA
44035
440-329-6000
800-333-6900
Technical Services
440-329-6593
800-832-4707

www.invacare.com



EU Representative

Invacare Deutschland GmbH
Kleiststraße 49
D-32457 Porta Westfalica
Germany
Tel: (49) (0) 5731 754 0
Fax: (49) (0) 5731 754 52191



Manufacturer

Invacare Corp.
2101 E. Lake Mary Blvd.
Sanford, FL 32773
Tel: 407-321-5630



1176148-D 2017-7-14



Making Life's Experiences Possible®



Yes, you can.®