
Ameda Pearl Electric Breast Pump
Tire-lait électrique Pearl d'Ameda
Bomba para pechos eléctrica Ameda Pearl

- EN Instructions for use | Page 4
- FR Mode d'emploi | Page 32
- ES Instrucciones de uso | Page 57

Congratulations on joining the Ameda family!

We are excited to help you meet your breastfeeding goals along your breastfeeding journey with your baby.

If you have any questions regarding your Ameda breast pump, our certified ParentCare team is available by phone during business hours on weekdays at **866.992.6332**.

For more information about breastfeeding, breast pumping and other Ameda products visit www.ameda.com. For more product information or feedback in other countries, visit www.ameda.com for a list of distributors outside of the U.S.

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


1.	WARNINGS & CAUTIONS	7
2.	INTENDED USE	9
3.	CONTRAINDICATIONS	9
4.	PRODUCT DESCRIPTION	10
5.	SET-UP INSTRUCTIONS	11
	5.1 Cleaning and Disinfecting the HygieniKit	11
	5.2 Assembling the HygieniKit Milk Collection System	12
6.	INSPECTING BEFORE PUMPING	13
7.	OPERATING INSTRUCTIONS	14
	7.1 Plug into Power Source	14
	7.2 Button Control	14
	7.3 LCD Icons & Modes	15
	7.4 To Begin Pumping	16
	7.5 How to Fit Your Flange	18
	7.6 Breast Milk Storage Guidelines	19
8.	AMEDA PEARL ELECTRIC BREAST PUMP TROLLEY	20
9.	CARRY CASE	21
10.	CLEANING THE PUMP	21
11.	MAINTENANCE	21
	12.1 Breast Pump Valve and Diaphragm Check	21
	12.2 Storing & Transporting	21
12.	TROUBLESHOOTING	22
13.	WARRANTY	23
14.	TECHNICAL DATA	24
	14.1 Specifications	24
	14.2 Electromagnetic Compatibility	25
	14.3 Product Disposal	29
15.	OPERATIONAL & STORAGE CONDITIONS	29
16.	TRAVEL	29
17.	AMEDA ACCESSORIES & SPARE PARTS	30
18.	SYMBOLS REFERENCE	31

Icon Glossary

	Warning		Wash hands		Serial number ⁹		Alternating current ¹⁴
	Caution		Temperature range ⁴		Catalogue number ¹⁰		Direct current ¹⁵
	Important		Humidity range ⁵		Unit of measure		On/Off button ¹⁶
	Follow instructions for use		Atmospheric pressure range ⁶		Class II electrical appliance (Double insulated) ¹¹		UL Safety Certification Mark
	Fragile, handle with care ¹		Expiration		Type BF applied part ¹²		For indoor use only ¹⁷
	Protect from sunlight ²		Manufacturer ⁷		Waste electrical and electronic equipment ¹³		Non-ionizing Electromagnetic Radiation ¹⁸
	Keep dry ³		Batch Code ⁸		Solid particle/liquid ingress classification		FCC Compliance

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

EN 1. WARNINGS & CAUTIONS

IMPORTANT SAFEGUARDS

When using electrical products, especially if children are present, the following basic safety precautions should always be followed.

The Ameda Pearl Electric Breast pump is for multiple-users in a hospital setting and for a single-user at home. Multiple-user pumps are designed so that the breast milk can never touch the working parts of the pump that are shared. The only part of a multiple-user breast pump that you can safely share is the pump unit.

For proper operation, the pump must be connected to the Ameda HygieniKit Milk Collection system breast pump kit (purchased or acquired separately). For your health and safety, the HygieniKit should never be shared between mothers or resold from one mother to another.

READ ALL INSTRUCTIONS BEFORE USING THIS PRODUCT

DANGER: To reduce the risk of electrocution:

- Always unplug electrical product immediately after use.
- **DO NOT** use while bathing, showering or swimming.
- **DO NOT** place or store where product can fall or be pulled into a bathtub, shower, sink, or pool.
- **DO NOT** place or drop product into water or other liquid.
- **DO NOT** reach for electrical product that has fallen into water. Immediately unplug from electrical outlet and contact manufacturer.
- **DO NOT** run water over breast pump.

The following symbols identify all instructions that are important to safety. Failure to observe these instructions can lead to injury to yourself or damage to the breast pump.

	WARNING Can lead to serious injury or death.
	CAUTION Can lead to minor injury.
	IMPORTANT Can lead to material damage.

WARNING: To avoid fire, electrocution, serious burns or injury:

- **NEVER** leave product unattended when plugged into an electrical outlet.
- Close supervision is necessary when product is used near children or persons with disabilities.
- Use product for its intended use **ONLY**, as described in this manual.
- **DO NOT** use parts that are not recommended by the manufacturer.
- Prior to use, always inspect power adapter. If damage and/or exposed wires are found, do not use power adapter. Call Ameda ParentCare at **1.866.992.6332**.
- **NEVER** operate electrical products with a damaged cord or plug, if not working properly, if dropped, damaged, or becomes wet.
- Keep cord away from heated surface.
- **NEVER** use while sleeping or when drowsy.
- **NEVER** insert or drop any object into any opening or into tubing.
- **KEEP AWAY** from heated surfaces or open flames, the breast pump and detachable components are not heat-resistant.
- **DO NOT** use near flammable materials.

Environmental Conditions That Affect Use

- **DO NOT** use outdoors with electrical cord.
- **DO NOT** operate where oxygen or aerosol spray products are being used.

WARNING: To avoid health risk and reduce the risk of injury:

Note: Do not use product if damaged. Use product only for intended use as described in this manual.

- **DO NOT** use pump while driving.
- Pumping can induce labor, **DO NOT** use while pregnant unless approved by your breastfeeding specialist or healthcare professional.
- **DO NOT** pump while sleeping or drowsy, it could result in tissue damage.
- After every use, clean all parts (flange, diaphragm, valve and bottle) that come in contact with your breast and breastmilk.
- **DO NOT** use attachments other than those recommended by manufacturer
- **DO NOT** make mechanical or electrical modifications to breast pump base.
- Before each use of the Ameda Pearl Electric Breast Pump and HygieneKit system, visually inspect the individual components for cracks, chips, tears, discoloration or deterioration. In the event that damage to the device or kit is observed, please discontinue use until the parts have been replaced. For the replaceable parts, please call your local distributor or the location where you purchased the pump or kit. For a listing of distributors in your country, please visit www.ameda.com.
- **DO NOT** repair or modify the device by yourself; contact Ameda if a problem occurs.
- **DO NOT** continue pumping for more than 2 consecutive pumping sessions if no results are achieved.
- **DO NOT** heat or thaw frozen breastmilk in a microwave or boiling water. Microwaving can heat liquids unevenly. Both microwaving and boiling water can result in possible burns.
- If breast pump kit is not sterile, clean and disinfect all parts that come into contact with your breast and breast milk prior to first use.
- If you are a mother who is infected with Hepatitis B, Hepatitis C or Human Immunodeficiency Virus (HIV), pumping breast milk will not reduce the risk of transmitting the virus to your baby through your breast milk.
- If irritation or discomfort occurs, discontinue use and see a doctor.
- Accessible materials may cause allergic reactions.
- **DO NOT** service the product while in use. The only user-replaceable components are the overlay panel and rubber feet.
- **DO NOT** add more than 5 lbs weight to trolley basket or assemble the trolley basket higher than 19 inches from the floor to avoid risk of tipping and falling of product.

STRANGULATION HAZARD: Keep tubing and power adapter cords out of reach of children or persons with disabilities.

CHOKING HAZARD: Small parts, **NOT** for children under 3 years.

CAUTION: Can lead to minor injury:

- **DO NOT** wrap cord around the body of the adapter.
- **USE** only the AC power adapter that is supplied with the Ameda Pearl Electric Breast Pump.
- **DO NOT** turn on the breast pump in an aircraft.
- **DO NOT** shorten tubing or modify the breast pump kit.
- Plug the AC power adapter into the pump base prior to plugging into the wall socket.
- Remove AC power adapter from power source before cleaning the pump.
- Keep away from dust, lint and direct sunlight to avoid potential reduction in pump performance.
- **DO NOT** attempt to remove the breast flange from your breast while pumping. If you are not feeling comfortable or feeling pain with pumping, turn the pump off, break the seal between the breast and breast flange with your finger, and remove breast flange from your breast gently.
- **ALWAYS** place pump on a horizontal flat surface when using product.
- **BATTERY CANNOT BE REMOVED.** This product contains a Lithium-ion battery. The battery or internal components cannot be replaced. **DO NOT** attempt to open pump base.
- Wash all parts that are exposed to breast milk immediately after use. This will help remove breast milk residue and prevent growth of bacteria.

IMPORTANT: Can lead to material damage.

- **NEVER** put electrical product in water or a sterilizer, as you can cause permanent damage to the breast pump.
- **DO NOT** use antibacterial or abrasive cleaners/detergents when cleaning breast pump or breast pump parts.
- **DO NOT** turn the breast pump on when the tubing port cover on the front of your pump is closed.
- When plastic bottles and plastic breast pump kit parts are frozen, they can become brittle and may break when dropped.
- **DO NOT** overfill bottles during pumping or when transferring to an empty bottle. Replace bottles when milk reaches top fill line.
- Handle carefully. Bottles and breast pump kit parts may become damaged if dropped, overtightened, or knocked over.
- Always check bottles and other breast pump kit parts for damage prior to pumping. Do not use those parts if damaged.

SAVE THESE INSTRUCTIONS

IMPORTANT NOTE for European Users: If a serious incident occurs in relation to the use of this breast pump, please report it immediately to Ameda at parentcare@ameda.com and to the competent authority of the Member State in which you are located. A list of vigilance contact points in the EU can be found at https://ec.europa.eu/health/sites/default/files/md_sector/docs/md_vigilance_contact_points.pdf.

2. INTENDED USE


The Ameda Pearl Electric Breast Pump is a powered breast pump to be used by lactating women to express and collect milk from their breasts. It is intended for multiple users in a hospital setting. It is also intended for home use by a single user.

3. CONTRAINDICATIONS

There are no known contraindications for use with this product.

4. PRODUCT DESCRIPTION

Ameda Pearl Electric Breast Pump & HygieniKit™ Milk Collection System

Description

1. Ameda Pearl Electric Breast Pump Unit

2. Power Adapter

Contents of an Ameda HygieniKit Milk Collection System

3. Tubing Adapter

4. Tubing

5. Adapter Cap

6. Diaphragm

7. Breast Flange (25.0mm)

8. Silicone Valve

9. Bottle

10. Bottle Locking Ring Cap

11. Bottle Locking Disc

The subject device is compatible with the HygieniKit Milk Collection System, but the HygieniKit Milk Collection System is not included in the package and is sold separately.

5. SET-UP INSTRUCTIONS

5.1 Cleaning and Disinfecting the HygieniKit

BEFORE EVERY USE

WASH HANDS THOROUGHLY with soap and water before touching breast pump, kit and breasts. Inspect and assemble a clean pump kit. If your tubing is moldy, discard and replace immediately.

WHEN USING A NEW NON-STERILE KIT

1. Inspect all kit components for damage (cracks, tears, etc.)
2. Place only the following parts in a pot of boiling water: diaphragm, flange, valve, bottles, bottle ring and locking discs
3. Boil for at least 10 minutes
4. Inspect all boiled parts again


WARNING: DO NOT leave boiling parts unattended. Boiled water and steam from boiled water can cause severe skin burns, so use caution when boiling parts in water and removing parts from boiling water. See the Center for Disease Control website <https://www.cdc.gov/healthywater/pdf/hygiene/breast-pump-fact-sheet.pdf> for additional cleaning guidelines.

AFTER EVERY USE


Store milk safely. Cap milk collection bottle or seal milk collection bag, label with date and time, and immediately place in a refrigerator, freezer, or cooler bag with ice packs.

Clean pumping area, especially if using a shared pump. Clean the pump as specified in this manual.

DO NOT wash adapter cap, tubing and tubing adapter.

Take apart breast pump tubing and separate all parts that come in contact with breast/breast milk.

Rinse breast pump parts that come into contact with breast/breast milk by holding under running water to remove remaining milk. **Do not place parts in sink to rinse.**

Clean pump parts that come into contact with breast/breast milk as soon as possible after pumping. You can clean your hard plastic parts (FLANGES AND BOTTLES) on the top rack of your dishwasher or by hand. ALWAYS CLEAN THE SOFT SILICONE PARTS (DIAPHRAGMS AND VALVES) BY HAND. When hand-washing any pump parts, always use a wash basin used only for cleaning the pump kit and infant feeding items.

Place pump parts in a clean wash basin used only for infant feeding items. **Do not place pump parts directly in the sink!**

Add soap and hot water to basin.


Scrub items using a clean brush used only for infant feeding items. Do not use a brush on valves or diaphragms.

Rinse by holding items under running water, or by submerging in fresh water in a separate basin. When rinsing the valves, use a GENTLE STREAM of running water. Valves may be damaged by the volume and pressure the faucet is capable of producing.

Air-dry thoroughly. Place pump parts, wash basin, and bottle brush on a clean surface in an area protected from dirt and dust. Do not use a dish towel to rub or pat items dry!

Clean wash basin and bottle brush. Rinse them well and allow them to air-dry after each use. Wash them by hand or in a dishwasher at least every few days.

For detailed instructions on cleaning or disinfecting your pump parts, visit <https://www.cdc.gov/healthywater/pdf/hygiene/breast-pump-fact-sheet.pdf>

IMPORTANT: The flanges and milk storage containers can also be washed on upper rack of a standard, household dishwasher. **DO NOT** use solvents or abrasives. Allow to air dry on a clean surface.

TUBING:

It is necessary to wash the tubing only if condensation is present. If condensation is present in the tube, wash and dry it immediately or replace the tubing. If milk is present in the tubing, do not attempt to wash or clean the tubing. Instead contact Ameda Customer Service at **1.866.992.6332**.


HOW TO WASH TUBING:

1. Rinse the tubing by pouring cool water and make sure the water flows through the entire tubing from one end to the other end.
2. Wash the tubing in warm, soapy water and rinse thoroughly.


HOW TO DRY TUBING: Hang the tubing to air dry and make sure that the tubing is completely dry before using.

5.2 Assembling the HygieniKit™ Milk Collection System


Once your pump parts are dry, wash your hands and assemble the HygieniKit Milk Collection System (Pump Kit).


HOLD by the rim or thicker sides of the valve.

PUSH the valve onto the bottom of the flange.


ATTACH and **TIGHTEN** the bottle or attach milk storage bag.

INSERT the diaphragm into the top of the flange.

It is important that the inside of the diaphragm is completely DRY to avoid moisture from forming in the tubing.


The tubing should still be **CONNECTED** to the adapter cap(s) and tubing adapter.


If you have taken these apart, attach tubing to adapter cap(s) and the tubing to the tubing adapter.

ATTACH the adapter cap(s) by pressing it onto the top of the flange(s) until you hear a pop.

ASSEMBLED HYGIENIKIT for single pumping.

ASSEMBLED HYGIENIKIT for double pumping.

6. INSPECTING BEFORE PUMPING


IMPORTANT: Before using the Ameda Pearl Electric Breast Pump and HygieniKit™ system, visually inspect the individual components for cracks, chips, tears, discoloration or deterioration. In the event that damage to the device or kit is observed, please discontinue use until the parts have been replaced. For the replaceable parts, please call your local distributor or the location where you purchased the pump or kit. For a listing of distributors in your country, please visit www.ameda.com.

7. OPERATING INSTRUCTIONS

7.1 Plug into Power Source (to charge battery or run on AC power)

Choosing and checking your electrical source

AC POWER ADAPTER


WARNING: Use **ONLY** the AC power adapter supplied with the Ameda Pearl Electric Breast Pump: Model No. UES36LCP1-150200SPA.

WARNING: If the power adapter housing or wiring becomes loose, separated, or frayed, stop use of the power adapter immediately and contact customer service at **1.866.992.6332**. Stop using the device immediately if you observe any smoke or burning from the pump unit or power adapter.

1. **PUSH** small end of the AC power adapter into the socket at the rear of the pump base as shown in section 7.4
2. **PLUG** other end into the AC outlet as shown above

7.2 Button Control

A. POWER ON/OFF BUTTON

B. MODE SELECTION BUTTON

Switches between Stimulation Mode and Expression Mode

C. - / + VACUUM LEVEL

Decreases or increases suction strength during pumping



D. - / + SPEED

Decreases or increases suction speed (cycles per minute)

E. BACKLIGHT BUTTON

F. NIGHTLIGHT BUTTON

7.3 LCD Icons & Modes

A. TIMER DISPLAY

Tracks pumping time. Automatically shuts-off after 60 minutes.

B. STIMULATION MODE

Mimics baby's fast suction pattern with fast speed and lower suction to start milk flow. This icon will disappear when in Expression Mode.

NOTE: After two minutes, the pump will automatically switch to Expression Mode unless you manually adjust beforehand.

C. SUCTION STRENGTH LEVEL

Level of suction during pumping. The Ameda Pearl Electric Breast Pump has two suction modes (Dual-Phase), each with adjustable suction levels:

- Stimulation Mode: 1 - 6 levels
- Expression Mode: 1 - 12 levels

Adjust the suction level by pressing the **-** or **+** button to decrease and increase level, respectively.

D. EXPRESSION MODE

Mimics baby's slow, steady suction pattern with slower speed and higher suction to express more milk gently and efficiently. This icon will disappear when in Stimulation Mode.

E. SPEED

Stimulation Mode: Fixed speeds at different suction levels as follows:

- Suction Strength Level 1 = 120 cycles per minute
- Suction Strength Level 2,3 = 100 cycles per minute
- Suction Strength Level 4, 5, 6 = 80 cycles per minute

Expression Mode: 6 individually selectable speeds at any Suction Strength Level

Adjust the speed level in expression mode by pressing the **-** or **+** button to decrease and increase level, respectively.

F. BATTERY LEVEL

Standard charge time, 0 to 100% capacity: At least 2 hours.

Battery operation time when the battery is fully charged: up to 4 hours.

Blinking battery segment = Battery charging.


Blinking battery status indicator = Low battery.

If the battery status indicator indicates one segment left, the pump will be off in 30 minutes.

7.4 To Begin Pumping


STEP 1: PUMP PLACEMENT

Place pump on a horizontal flat surface.


STEP 2: CONNECT THE POWER ADAPTER

Insert the power adapter into the rear of the pump and an AC power outlet. If operating on battery alone, use of the power adapter is not required.

STEP 3: INSERT TUBING


If the tubing port cover is not open, grasp the cover near the top and gently pull the cover open and expose the tubing port. Insert the tubing adapter into the front of pump base, then turn slightly to right while pushing firmly.


Single v. Double Pumping

DOUBLE PUMPING = expressing from both breasts at the same time


SINGLE PUMPING = expressing from just one breast at a time

A. For DOUBLE PUMPING:
Both tubes should be attached to the tubing adapter.

B. For SINGLE PUMPING:
Remove one tube from the tubing adapter and close one side using the attached plug.


STEP 4: PLACE THE FLANGE ON BREAST

Center flange on breast and press firmly to create a tight seal.


Check breast flange fit, if pumping causes discomfort a larger or smaller breast flange size may be needed. See section on how to properly fit your flange. See section 7.5 below on how to properly fit your flange.

NOTE: Make sure your upper body is upright.

STEP 5: POWER ON THE PUMP

Press the Power button to turn on the device. It will begin in Stimulation Mode.


STEP 6: ADJUST STIMULATION VACUUM LEVEL

Adjust the Stimulation Vacuum Level between 1 and 6 to your preferred comfort level. Never set to an uncomfortable level.


TIP: Start off with the suction low, and then gradually increase the suction setting to your comfort level. You may find that the level of comfortable suction varies from day to day, during the day and even during the same pumping session.

NOTE: After two minutes, the pump will automatically switch to Expression Mode unless adjusted manually.

STEP 7: CHANGE MODE


Once your milk starts to flow using the STIMULATION MODE, press the Mode button to switch into Expression Mode. Once your milk stops flowing freely while in Expression Mode, press the Mode button to switch back to Stimulation Mode to induce a second let-down if desired.

STEP 8: ADJUST EXPRESSION VACUUM LEVEL


Adjust the Expression Vacuum Level between 1 and 12 to your preferred comfort level. Never set to an uncomfortable level.

TIP: Start off with the suction low, and then gradually increase the suction setting to your comfort level. You may find that the level of comfortable suction varies from day to day, during the day and even during the same expression session.

STEP 9: ADJUST SPEED


Adjust the speed between 1 and 6 to your preferred comfort level. Never set to an uncomfortable level.

STEP 10: END PUMPING SESSION


When you wish to stop expressing milk, first press the Power button, then remove milk collection kit from your breast.

TIP: After pump is off, if necessary, INSERT a finger between breast and breast flange to break suction.

SECURE BOTTLES


Place bottles in center of bottle holder compartment(s) in pump.

If disconnecting the milk collection kit from the pump, close the tubing port cover by pressing cover firmly over the tubing port.


7.5 How to Fit Your Flange

Check breast flange fit. If pumping causes discomfort, a larger or smaller breast flange size may be needed.


GOOD FIT

During pumping, your nipple moves freely in the breast flange tunnel. You see space around the nipple. Not much areola is drawn into the tunnel with the nipple.

FLANGE TOO SMALL

During pumping, some or all of your nipple rubs against the sides of the breast flange tunnel.

FLANGE TOO LARGE

During pumping, more areola is drawn into the breast flange with your nipple. Your areola may rub against the side of the breast flange tunnel.

7.6 Breast Milk Storage Guidelines

- 1: Write the date on the storage container. Include your child's name if you are giving the milk to a child care provider.
- 2: Gently swirl the container to mix the cream part of the breastmilk that may rise to the top back into the rest of the milk. Do not shake the milk as this can make some of the milk's valuable parts break down.
- 3: Refrigerate or chill milk right after it is pumped, if possible. You can put it in the refrigerator, place it in a cooler or insulated cooler pack, or freeze it in small (2- to 4-ounce) batches for later feedings. Pumped milk is OK without being refrigerated for up to 4 hours after pumping.

	Room Temperature 77°F (25°C)	Time in Refrigerator 39°F (4°C)	Time in Freezer 0°F (-18°C)
Freshly Expressed or Pumped	Up to 4 Hours	Up to 4 Days	6 Months - Best 12 Months - Acceptable
*Thawed, Previously Frozen	1-2 Hours	Up to 1 Day (24 Hours)	Do Not Refreeze
Leftover From Feeding (baby did not finish bottle)	Use within 2 hours after the baby is finished feeding		

***NOTE:** When thawing frozen milk, label as thawed when completely thawed (i.e. no ice crystals present). Use the time when completely thawed to base acceptable time limits for use rather than when it is taken from the freezer. These guidelines for storage and thawing of breast milk are a recommendation. Contact your lactation consultant or breastfeeding specialist for further information.

These guidelines are for healthy full-term babies and may vary for premature or sick babies. Check with your health care provider. Find more breastfeeding resources at: WICBreastfeeding.fns.usda.gov and www.cdc.gov/breastfeeding.

SAFE THAWING OF BREAST MILK

Always thaw the oldest breast milk first. Remember first in, first out. Over time, the quality of breast milk can decrease.

There are several ways to thaw your breast milk:

- In the refrigerator overnight.
- Set in a container of warm or lukewarm water.
- Under lukewarm running water.
- Never thaw or heat breast milk in a microwave. Microwaving can destroy nutrients in breast milk and create hot spots, which can burn a baby's mouth.
- If you thaw breast milk in the refrigerator, use it within 24 hours. Start counting the 24 hours when the breast milk is completely thawed, not from the time when you took it out of the freezer.
- Once breast milk is brought to room temperature or warmed, use it within 2 hours.
- Never refreeze breast milk after it has thawed.
- Additional information: https://www.cdc.gov/breastfeeding/recommendations/handling_breastmilk.htm

8. AMEDA PEARL ELECTRIC BREAST PUMP TROLLEY

The Ameda Pearl Breast Pump Trolley is not included with the Ameda Pearl Electric Breast Pump and is sold separately.

Attaching the Ameda Pearl Electric Breast Pump to Trolley (Trolley is an optional accessory)

NOTE: Attach pump to trolley prior to plugging cord into wall outlet.

- 1: When assembling Trolley, ensure not to place the basket any higher than 19 inches from the floor (height to bottom of basket).
- 2: Place the three rubber washers on the trolley mounting plate to align over the three holes (C).
- 3: Place the Ameda Pearl Breast Pump on the trolley mounting plate so that the basket lies underneath the control panel. Ensure the rubber feet (A) are on the mounting plate and the three holes (B) in the bottom of the pump are aligned with the three washers and holes (C) in the mounting plate.
- 4: Use three mounting screws supplied with the trolley to secure the pump to the trolley mounting plate.


A. Rubber Feet (4)

B. Pump Mounting Holes (3)

C. Trolley Plate Mounting Holes (3)


D. Trolley Basket (1)

Pump Bottom View:

**Breast Pump
Bottom Front**

Trolley Top View:

Trolley Mounting Plate Front

Trolley Operating Instructions

- 1: Once pump is attached to trolley, ensure all locking wheels are unlocked to move the trolley-pump assembly. Wheels are unlocked by lifting the locking tab on each locking plate.
- 2: Move the trolley-pump assembly by gently pushing.
- 3: When trolley-pump assembly is in the desired position, lock all locking wheels by pushing down on the locking tab on each locking wheel.

9. CARRY CASE

Carry Case is not included with the Ameda Pearl Electric Breast Pump and is sold separately.

The Carry Case is designed to hold (1) Ameda Pearl Electric Breast Pump, (1) AC power adapter, and (1) HygieniKit breast milk collection kit.

Place the Ameda Pearl Electric Breast Pump and power adapter into the spaces designed for them in the Carry Case to reduce movement while transporting the pump.

To close case push latches down until latches snap into place. Confirm both latches are completely closed and secured before carrying case. To open case lift latches up.

10. CLEANING THE PUMP

Turn the breast pump off and unplug the pump from the power source.

For at-home/end-users

- After each pump session, wipe the Ameda Pearl Electric Breast Pump with a clean, damp cloth.


For Hospitals/Rental Stations

- Clean between users.
- Wear disposable gloves and use a hospital approved antibacterial cleaner such as CaviCide®, CaviWipes®, or PDI® Sani-Cloth® AF3.
- Wipe Ameda Pearl Electric Breast Pump clean and remove all surface debris from the device, including the power cord, trolley (if used) and carry case.
- Repeat cleaning step as necessary using a new cloth or wipe each time.
- When finished cleaning, allow cleaned surfaces to air dry for a minimum of 10 minutes.

11. MAINTENANCE

The breast pump base does not require service or maintenance. The user may replace the overlay panel and rubber feet as needed, using Ameda supplied components.

11.1 Breast Pump Valve and Diaphragm Check

Inspect the valves and diaphragms for tears or holes. Replace any damaged valves and diaphragms prior to pumping. Wash valves and diaphragms per section 8 "Cleaning the HygieniKit". Do not insert bottle brush or something similar into valves. We recommend replacing valves every 8-12 weeks and diaphragms every 4-6 weeks. When pumping at maximum settings, you may need to replace diaphragms more frequently.

11.2 Storing & Transporting

Disconnect HygieniKit and power adapter from pump before transporting or storing the product. Be sure to close the tubing port cover by pressing cover firmly over the tubing port. Handle with care when transporting product. Store product in cool, dry place and protect from dust and dirt. Breast Pump and power adapter may also be stored inside the Ameda Pearl Electric Breast Pump carry case.

12. TROUBLESHOOTING

<p>SLOW MILK FLOW, UNEVEN OR NO SUCTION</p>	<p>VALVE AND DIAPHRAGM CHECK</p> <p>Inspect the valves and diaphragms for tears or holes.</p> <p>Replace any damaged valves and diaphragms prior to pumping. Wash valves and diaphragms per cleaning section. Do not insert bottle brush or similar into valves. We recommend replacing valves every 8-12 weeks and diaphragms every 4-6 weeks but this can vary depending on cleaning and use.</p> <p>CONNECTION CHECK</p> <p>Check assembly is correct and connection points are securely attached.</p> <ul style="list-style-type: none"> • Tubing adapter to tubing port • Tubing adapter to tubing <p>Check adapter cap is attached to breast flange. It may be useful to remove adapter cap and reattach to breast flange.</p> <p>FLANGE/BOTTLE CHECK</p> <p>Check flanges, adapter cap or bottles for chips, cracks, breaks, or other damage. Replace any damaged flanges/adapter caps/bottles prior to pumping.</p> <p>BREAST FLANGE FIT CHECK</p> <ul style="list-style-type: none"> • Good Fit: Nipple is moving freely in tunnel of flange. • Too Small: Nipple is rubbing against tunnel of flange. • Too Large: Areola is being pulled into tunnel of flange. • Check to make sure the breast shields form a complete seal around the breast while pumping.
<p>AUDIBLE, HISSING NOISE OR AIR LEAKAGE</p>	<p>VALVE AND DIAPHRAGM CHECK</p> <p>Inspect the valves and diaphragms for tears or holes.</p> <p>Replace any damaged valves and diaphragms prior to pumping. Wash valves and diaphragms per cleaning section. Do not insert bottle brush or similar into valves. We recommend replacing valves every 8-12 weeks and diaphragms every 4-6 weeks but this can vary depending on cleaning and use.</p> <p>CONNECTION CHECK</p> <p>Check assembly is correct and connection points are securely attached.</p> <ul style="list-style-type: none"> • Tubing adapter to tubing port • Tubing adapter to tubing <p>Check adapter cap is attached to breast flange. It may be useful to remove adapter cap and reattach to breast flange.</p>
<p>POWER SOURCE</p>	<p>AC POWER ADAPTER</p> <p>Make sure the AC power adapter is plugged all the way into the rear of your pump base. Firmly press the power button on the lower left side of the pump face. Make sure the outlet is working.</p>
<p>"Lo Batt" Notification on LCD Screen</p>	<p>This indicates that the battery level is low. Use the AC Power Adapter to recharge the battery and operate the pump.</p>
<p>"Check Kit" Notification on LCD Screen</p>	<ul style="list-style-type: none"> • Ensure tubing adapter is properly fitted into tubing adapter port on pump. • Check that adapter cap is firmly snapped onto top of breast flange. • Ensure silicone diaphragm is in place inside top of breast flange. • Check tubing and silicone diaphragm to ensure that parts are free of defects. • Check tubing is firmly attached to adapter cap. • If problem persists, contact Ameda at 1.866.992.6332 • Check to make sure the breast shields form a complete seal around the breast while pumping.

If there are still questions or issues, please call Ameda's certified customer service ParentCare team at **1.866.992.6332**.

13. WARRANTY

THE AMEDA WARRANTY IS NONTRANSFERABLE BETWEEN INDIVIDUALS. THE AMEDA PEARL ELECTRIC BREAST PUMP ("PRODUCT") IS WARRANTED TO THE ORIGINAL USER – WHETHER INSTITUTIONAL OR INDIVIDUAL ("FIRST USER") – ONLY.

From the date of the Individual User's purchase of this product, Ameda, Inc. ("Ameda") warrants the Product to the Individual User against defects in material or workmanship for a period of three years for the pump mechanism and 90 days on the HygieniKit except for breast pump diaphragms and valves. Ameda's sole obligation under this limited warranty shall be to replace, at Ameda's option, any Product that is determined to be defective by Ameda and determined to be covered by this express limited warranty. Replacement under this limited warranty is the sole and exclusive remedy of the First User.

PROOF OF PURCHASE IN THE FORM OF A RECEIPTED INVOICE OR BILL OF SALE EVIDENCING THAT THE PRODUCT IS WITHIN THE WARRANTY PERIOD MUST BE PRESENTED TO OBTAIN WARRANTY SERVICE. This limited warranty is extended by Ameda ONLY to the First User and is not assignable or transferable. For warranty assistance contact Ameda at **1.866.992.6332**.

EXCEPT TO THE EXTENT PROHIBITED BY APPLICABLE LAW, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE ON THE PRODUCT IS HEREBY DISCLAIMED. AMEDA SHALL NOT BE LIABLE FOR ANY INCIDENTAL, CONSEQUENTIAL, SPECIAL OR PUNITIVE DAMAGES OF ANY KIND FOR BREACH OF THE EXPRESS LIMITED WARRANTY ON THE PRODUCT OR ANY WARRANTY IMPLIED BY OPERATION OF LAW. OTHER THAN THE EXPRESS LIMITED WARRANTY SET FORTH ABOVE, THERE ARE NO OTHER WARRANTIES THAT ACCOMPANY THIS PRODUCT AND ANY ORAL, WRITTEN OR ANY OTHER REPRESENTATION OF ANY KIND IS HEREBY DISCLAIMED.

NON-APPLICABILITY OF WARRANTY

In addition to any other limitations on the warranty discussed above, the warranty shall expire and be of no force or effect upon the occurrence of any of the following:

- (i) the cause of damage to the Product is due to improper handling or use.
- (ii) opening of the pump motor, modifications, self-repair or use outside the instructions in this manual voids warranty.
- (iii) the Product is not properly maintained by the Individual User as determined by Ameda.
- (iiii) non-Ameda approved accessories or spare parts are used with the Product as determined by Ameda.

If any of the above occurs, Ameda shall not be responsible for damage to the Product, parts thereof or injuries arising therefrom, either directly or indirectly. The warranty terms set forth above replaces all prior warranties with respect to the purchased Product, whether in writing or otherwise.

Outside the United States: Call your local distributor or location where you purchased the Product. For a listing of distributors in your country, please visit www.ameda.com.

14. TECHNICAL DATA

14.1 Specifications

Product Description	Electric Breast Pump
Model	Ameda Pearl
Suction	Stimulation mode: Approx. 50 to 150mmHg
	Expression mode: Approx. 30 to 250mmHg
Vacuum Levels	Stimulation mode: 6 Levels
	Expression mode: 12 Levels
Cycles Per Minute	Stimulation mode: 80 to 120 CPM
	Expression mode: 30 to 48 CPM
Power Source	Medical AC Adapter - Output: 15VDC, 2A, 30W. Input: 100-240VAC, 50/60Hz, 1A.
	2600mAh Lithium-ion rechargeable battery
Battery Life	At least 3 hours after full charge cycle
Timer Function	Automatic Power-off after 60 minutes
Unit Weight	Approx. 1.52 kg (3.35 lbs)
Unit Dimensions	110 x 93 x 55mm (L x W x H)
Operating Environment	Temperature: 5°C to 40°C (41°F to 104°F)
	Humidity: 15% to 90%RH
	Pressure: 90kPa to 106kPa
Storage Environment Transportation	Temperature: -20°C to 55°C (-4°F to 131°F)
	Humidity: 10% to 90%
Expected Service Life	3 years
Ingress Protection Rating	IP22
Classification	Internal Powered Equipment, Type BF Applied Part
Protection Class	Class II/Double Insulated
Safety Compliance	UL CERTIFIED E49355 MEDICAL - GENERAL MEDICAL EQUIPMENT AS TO ELECTRICAL SHOCK, FIRE AND MECHANICAL HAZARDS ONLY IN ACCORDANCE WITH ANSI/AAMI ES 60601-1 (2005) + AMDI (2012) AND CAN/CSA-C222 NO. 60601-1 (2008) + (2014) AND IEC 60601-1-6 (2010) + AMDI (2013) AND IEC 62366 (2007) + AMDI (2014) AND IEC 60601-1-11:2015
	
	Device complies with the requirements of the Federal Communications Commission

Specifications are subject to change without notice.

14.2 Electromagnetic Compatibility

IMPORTANT INFORMATION REGARDING ELECTROMAGNETIC COMPATIBILITY (EMC) This electrical medical equipment needs special precautions regarding EMC and shall be put into service according to the EMC information provided in this user manual. The equipment conforms to IEC 60601-1-2 for both immunity and emissions. Nevertheless, special precautions need to be observed. This equipment is intended for Home Healthcare environment.

The essential performance of the Ameda Pearl Electric Breast Pump is to maintain user-selected levels of pump speed and vacuum pressure not to exceed 250mmHg + 25mmHg. The pump was tested for immunity to electromagnetic disturbances and passed according to this criteria.

WARNING: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

WARNING: Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

TABLE 1

Guidance and declaration of manufacturer-electromagnetic emissions		
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
Radiated emissions CISPR 11	Group 1, Class B	The device uses RF energy only for its internal function. Therefore, its emissions are very low and are not likely to cause any interference in nearby electronic equipment.
Conducted emissions CISPR 11	Group 1, Class B	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

FCC STATEMENT

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Caution: changes or modifications to this device not expressly approved by Ameda may void the user's authority to use the device.

Note:

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Note:

This equipment has been tested and found to comply with the immunity test levels in AIM 7351731 Rev. 2.0 (2017), Medical Electrical Equipment & System Electromagnetic Immunity Test for RFID Readers.

TABLE 2

Guidance and declaration of manufacturer-electromagnetic immunity (For home healthcare environment)			
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment			
IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment -guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transients/bursts IEC 61000-4-4	±2 kV, 100kHz, for AC power port	±2 kV, 100kHz, for AC power port	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±0.5kV, ±1kV (differential mode)	±0.5kV, ±1kV (differential mode)	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% UT; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% UT; 1 cycle and 70% UT; 25/30 cycles Single phase: at 0° 0% UT; 250/300 cycle	0% UT; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% UT; 1 cycle and 70% UT; 25/30 cycles Single phase: at 0° 0% UT; 250/300 cycle	Mains power quality should be that of a typical commercial or hospital environment.
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	30 A/m; 50Hz or 60Hz	30 A/m; 50Hz or 60Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

TABLE 3

Guidance and declaration of manufacturer-electromagnetic immunity (For home healthcare environment)			
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment			
IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3V for 0.15 - 80MHz; 6V in ISM and amateur radio bands between 0.15 - 80MHz 80% AM at 1kHz	3V for 0.15 - 80MHz; 6V in ISM and amateur radio bands between 0.15 - 80MHz 80% AM at 1kHz	Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Ameda Pearl, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.
Radiated RF IEC 61000-4-3	10V/m, 80MHz - 2.7GHz, 80% AM to 1kHz	10V/m, 80MHz - 2.7GHz, 80% AM to 1kHz	Interference may occur in the vicinity of equipment marked with the following symbol:
Proximity fields from RF wireless communications equipment IEC 61000-4-3	According to Table 9 of IEC 60601-1-2:2020	According to Table 9 of IEC 60601-1-2:2020	
Proximity magnetic fields IEC 61000-4-39	According to Table 11 of IEC 60601-1-2:2020	According to Table 11 of IEC 60601-1-2:2020	

TABLE 4

IEC 60601-1-2:2020 Table 11 - Test specifications for ENCLOSURE PORT IMMUNITY to proximity magnetic fields

Test frequency	Modulation	Immunity Test Level (A/m)
30 kHz ^{a)}	CW	8
134,2 kHz	Pulse modulation ^{b)} 2,1 kHz	65 ^{c)}
13,56 MHz	Pulse modulation ^{b)} 50 kHz	7,5 ^{c)}

a) This test is applicable only to ME EQUIPMENT and ME SYSTEMS intended for use in the HOME HEALTHCARE ENVIRONMENT.

b) The carrier shall be modulated using a 50 % duty cycle square wave signal.

c) r.m.s., before modulation is applied.

TABLE 5

IEC 60601-1-2:2020 Table 9 - Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment

Test frequency (MHz)	Band ^{a)} (MHz)	Service ^{b)}	Modulation ^{c)}	IMMUNITY TEST LEVEL (V/m)
385	380-390	TETRA 400	Pulse modulation ^{b)} 18 Hz	27
450	430 - 470	GMRS 460, FRS 460	FM ^{c)} ±5 kHz deviation 1 kHz sine	28
710	704 - 787	LTE Band 13, 17	Pulse modulation ^{b)} 217 Hz	9
745				
780				
810	800 - 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation ^{b)} 18 Hz	28
870				
930				
1720	1700 - 1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation ^{b)} 217 Hz	28
1845				
1970				
2450	2400 - 2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation ^{b)} 217 Hz	28
5240	5100 - 5800	WLAN 802.11 a/n	Pulse modulation ^{b)} 217 Hz	9
5500				
5785				
<p>NOTE If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the ME EQUIPMENT OR ME SYSTEM may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3.</p>				
<p>a) For some services, only the uplink frequencies are included. b) The carrier shall be modulated using a 50% duty cycle square wave signal. c) As an alternative to FM modulation, the carrier may be pulse modulated using a 50% duty cycle square wave signal at 18 Hz. While it does not represent actual modulation, it would be worst case.</p>				

WARNING: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Ameda Pearl Electric Breast Pump, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

WARNING: If the use location is near (e.g. less than 1.5 km from) AM, FM or TV broadcast antennas, before using this equipment, it should be observed to verify that it is operating normally to assure that the equipment remains safe with regard to electromagnetic disturbances throughout the expected service life. When the AC input voltage is interrupted, the equipment will stop battery charging and if the power supply restored, it could be recovered automatically. This degradation could be accepted because it will not lead to unacceptable risks and it will not result in the loss of basic safety or essential performance.

14.3 Product Disposal

1. The internal components of this product may contain hazardous materials. Dispose of this product according to local or regional waste administration systems and regulations. Caution: Hazardous materials present in electric and electronic equipment may harm human health and the environment if the equipment is not properly disposed of.
2. Do not dispose of electric or electronic equipment together with unsorted municipal waste.
3. For further information regarding the environmental performance of this product please visit our website www.ameda.com. The breast pump and parts are made of plastics which have an impact on the environment when disposed of. Please recycle according to local regulations.

15. OPERATIONAL & STORAGE CONDITIONS

Operational Conditions

Transport/Storage Conditions

16. TRAVEL

WARNING: Use **ONLY** the AC power adapter supplied with the Ameda Pearl Electric Breast Pump.

WARNING: Traveling by air with your pump. The lithium-ion battery contained in the Ameda Pearl Electric Breast Pump has a nominal watt-hour rating of 28 Wh. If you travel by air and plan to take the pump with you, be sure to check airline policy for any restrictions. If the airline allows the pump to be carried in checked baggage, make sure that the pump is completely switched off and protected against unintentional activation or damage. Use of the original packaging or the Ameda Pearl Electric Breast Pump carrying case is advised.

If you will be traveling outside of the country where you received your Ameda Pearl Electric Breast Pump, you will need to use an Ameda Pearl AC prong adapter specific to your destination in order to use your breast pump while traveling.

17. AMEDA ACCESSORIES AND SPARE PARTS

The AC power adapter is included with the purchase of the Ameda Pearl Electric Breast Pump. All other products shown below are sold separately. You may also be able to order these additional accessories or spare parts online, at a local Ameda retailer or Ameda insurance covered provider, visit www.ameda.com for more information.

CleanEase
Microwavable
Steam Sanitizing
Bag

AC Power
Adapter

One-Hand
Breast Pump

Bottle
Holders

21.0mm Flexishield™
Areola Stimulator

Breast Milk
Storage Bottles

Pump'N Protect
Bags

MiniFil Colostrum
Collection
Container 50mL

MiniFil Colostrum
Collection
Container 80mL

25.0mm
Breast Flange

28.5mm
Breast Flange

30.5mm
Breast Flange

Valve

Diaphragm

Tubing

Tubing
Adapter

Adapter
Cap

Locking Ring

Locking
Disk

Carry Case

Trolley

Cool'N Carry Milk
Storage Tote

Cooling
Element

18. SYMBOLS REFERENCE

1. ISO 15223-1, Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied, Part 1: General requirements, Clause 5.3.1, Fragile, handle with care / ISO 7000-0621, Graphical symbols for use on equipment, Fragile, handle with care
2. ISO 15223-1, Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied, Part 1: General requirements, Clause 5.3.2 Keep away from sunlight / ISO 7000-0624, Graphical symbols for use on equipment, Keep away from sunlight
3. ISO 15223-1, Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied, Part 1: General requirements, Clause 5.3.4, Keep dry / ISO 7000-0626, Graphical symbols for use on equipment, Keep away from rain
4. ISO 15223-1, Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied, Part 1: General requirements, Clause 5.3.7 Temperature Limit / ISO 7000-0632, Graphical symbols for use on equipment, Temperature Limit
5. ISO 15223-1, Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied, Part 1: General requirements, Clause 5.3.8 Humidity limitation / ISO 7000-2620, Graphical symbols for use on equipment, Humidity Limit
6. ISO 15223-1, Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied, Part 1: General requirements, Clause 5.3.9 Atmospheric pressure limitation / ISO 7000-2620, Graphical symbols for use on equipment, pressure Limit
7. ISO 15223-1, Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied, Part 1: General requirements, Clause 5.1.1 Manufacturer
8. ISO 15223-1, Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied, Part 1: General requirements, Clause 5.1.5 Batch code / ISO 7000- 2493, Graphical symbols for use on equipment, batch code
9. ISO 15223-1, Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied, Part 1: General requirements, Clause 5.1.7 Serial number / ISO 7000-2498, Graphical symbols for use on equipment, Serial number
10. ISO 15223-1, Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied, Part 1: General requirements, Clause 5.1.6 Catalogue number / ISO 7000- 2493, Graphical symbols for use on equipment, Catalogue number
11. IEC 60601-1, Medical electrical equipment – Part 1: General Requirements for basic safety and essential performance, Table D.1 Symbol 9 Class II equipment
12. IEC 60601-1, Medical electrical equipment – Part 1: General Requirements for basic safety and essential performance, Table D.1 Symbol 20 Type BF applied parts
13. EN 50419, Marking of Electrical and Electronic Equipment in accordance with Article 11(2) of Directive 2002/96/EC (WEEE)
14. IEC 60601-1, Medical electrical equipment – Part 1: General Requirements for basic safety and essential performance, Table D.1 Symbol 4 Alternating current
15. IEC 60601-1, Medical electrical equipment – Part 1: General Requirements for basic safety and essential performance, Table D.1 Symbol 1 Direct current
16. IEC 60601-1, Medical electrical equipment – Part 1: General Requirements for basic safety and essential performance, Table D.1 Symbol 29 Stand-by
17. IEC 60417-5957, Graphical symbols for use on equipment, Indoor Use only
18. IEC TR 60878:2015, Graphical symbols for electrical equipment in medical practice, Non-ionizing electromagnetic radiation

To register your pump: Please visit ameda.com/register-your-pump and input the serial number (SN) from the bottom of your pump.

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To register your pump: Please visit www.ameda.com and input the serial number (SN) from the bottom of your pump.

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Mothers Choice Products Ltd.
2133-11871 Horseshoe Way
Richmond, BC V7A5H5
1.800.604.6225

Distributed in Australia by:
Distribué en Australie par :
Distribuido en Australia por:
Midmed Pty Ltd.
Unit 4, 62 Borthwick Avenue
Murarrie, QLD 4172
61 (7) 3348 9155

Distribution in other countries: For product information or feedback, call your local distributor or location where you purchased the product. For a listing of distributors in your country, please visit www.ameda.com.

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Ameda Inc.